


JointMedica

Polymotion® Hip Resurfacing (PHR®): Size Chart

When carrying out a Polymotion® Hip Resurfacing operation, please make sure that the Acetabular Cup and Femoral Head are compatible.

You must use a Cup and a Head with the same corresponding letter.

Please see below for compatible sizes.



Letter for Compatibility	Femoral Head Size mm	Thin Wall Acetabular Cup 3mm	Thick Wall Acetabular Cup 4mm
A	40	46	48
B	42	48	50
C	44	50	52
D	46	52	54
E	48	54	56
F	50	56	58
G	52	58	60
H	54	60	62
I	56	62	64
J	58	64	66

Manufactured by

JointMedica

Distributed by

exactech

JointMedica Ltd, Unit 3 Ball Mill Top, Main Road, Hallow, Worcestershire, WR2 6LS. United Kingdom.

Tel: +44 (0)1905 640008 E-mail: cs@jointmedica.com